

Entrix 7

Low Profile Alternating Mattress Replacement

An alternating air mattress replacement system that offers a high degree of comfort and value for money.

Entrix 7 Mattress Replacement

- Sixteen 5-inch tall air cells.
- Quick release CPR rapidly deflates the entire mattress.
- Micro-ventilated air cells for great air circulation, optimal patient temperature regulation from shoulder to waist.
- Two way stretch coverlet washable, fully water proof and excellent moisture permeability.
- Built in 2 inches of foam no need for a separate mattress, easier for care staff to make beds, overall mattress height is lower.

Entrix Power Unit

- Dynamic/ Static modes.
- · Adjustable weight setting.
- IPS an integrated electronic pressure sensor for better user comfort, lower interface pressure and excellent results.
- Visual and audible alarms for power failure, low pressure.
- · Digital easy to use/ easy to clean.





Entrix 7 Product Specification

Individual Model No.

Power control Unit - S2202-3

Mattress: SR2203-2 - Standard

SR2203-2 - Narrow

Power Input REF SR366 Entrix

AC100 - 120V/ 50 & 60 Hz

AC220 - 240V/50 Hz

Power Consumption

Max. 10W

Power Unit Dimension

- 29.5 x 22.5 x 12.0 cm

Mattress Dimension

King Single - 200 x 107 x 18 cm (inflated)

- 200 x 90 x 18 cm (inflated)

Narrow

- 200 x 82 x 18 cm (inflated)

Material

Power Unit - ABS fire retardant UL Mattress Cell - 100% nylon with TPU

lamination

Coverlet - 100% nylon with PU backing

Maximum Weight capacity

150kg

Weight

Power control Unit - 2.5kg Mattress Set - 5kg

Product Compliance

EN60601 - 1 EN60601 - 2 UL60601 -1 HCPCS E0277

Operating & Transportation & **Storage Conditions**

Temperature - 5°C-45°C Humidity - 15%-60%

Entrix 7 & King Single					
Application	Risk	Weight Capacity	Mattress Type	TGA Registered	Warranty
Nursing Home, Home	Medium	30-150 kg (min-max)	Replacement	YES	2 YEARS



Micro ventilated air cells: 5" individual replaceable air cells with air ventalation from shoulder to waist



Quick release coupling; Easy connection between power unit and mattress for fast patient transferring



Fast deflation CPR facility

For pressure mapping information or to find a distributor in your state simply call 1300 732 004 or email natalie.nati@carilexmedical.com. For more information please visit www.carilexmedical.com.

Carilex Medical ANZ PO 613 Cherrybrook NSW 2126 P 1300 732 004 F 02 8834 0764 E natalie.nati@carilexmedical.com